

Chem 101 Laboratory Exercise #1 Laboratory Notebook

Measurement of Volumes and Weight: Accuracy and Precision

Using Microsoft Word, students are to write their in-lab notes below while completing the Laboratory exercise. The Laboratory notebook must be uploaded to the CHEM 101 Lab Brightspace site as a readable .pdf by the end of the lab period. Please see page 7 of the 202309 CHEM 101 lab manual as to all the information required in the in-lab notes. The documentation below must reflect the student's own work without assistance from others.

Name: **Arfaz Hossain** Lab Section: **B12** Quad: **2** Date: **21 September 2023**

In-lab Notes:

Experimental Procedures:

- 1 For weighing accurately an approximate 1.5 g of NaCl, I used a weigh boat to take the NaCl and then put that (subtracting the weight of the weigh boat) on a balance. The weight is documented in the Final Documentation.
- 2 Then I transferred the NaCl to a 50.00 mL volumetric flask from the weigh boat and filled to the mark with distilled water.
- 3 I took a clean dry 50 mL beaker on a balance, documented its weight (1.1).
- 4 Then I used a 10.00 mL volumetric pipette and transferred 10.00 mL of the NaCl solution from the volumetric flask to the beaker.
- 5 The weight of the solution without the beaker can be found in the Final Documentation (1.2).

I re-did the experiment for 4 more times.

Final Measurements:

Weight of the NaCl: **1.543g**.

1.1 Weight of the Beaker is **28.501g**.

1.2 Weight of the Solution-01: **11.166g**.

2.1 Weight of the Beaker is **28.503g**.

2.2 Weight of the Solution-02: **10.156g**.

3.1 Weight of the Beaker is **28.506g**.

3.2 Weight of the Solution-03: **10.183g**.

4.1 Weight of the Beaker is **28.508g**.

4.2 Weight of the Solution-01: **10.151g**.

Abstract:

By measuring the volume and weight of samples the density of a solution of NaCl was determined it to be **1.543**.

Data/Results:

10.00 mL of NaCl Solution	#1	#2	#3
Weight of sample from Volumetric Pipette			
Calculated Density			
Average Calculated Density of NaCl Solution			
%RSD for the density of NaCl Solution			

Algebraic Equations:

Density (g/mL) = weight of NaCl solution /volume of NaCl solution

Discussion:

The density of the prepared NaCl solution was determined to be ___ by measuring 3 samples of 10.00mL of the NaCl solution. The %RSD calculated for these measurements was ___. The accuracy of the measurement required the use of a pipette and not a cylinder because ___.

Conclusions:

The measured density of the NaCl solution using a volumetric pipette is ___ with a % relative standard deviation of ___.

References:

1. Reimer, M. et al, Laboratory Manual, Chemistry 101, pp. 13-18. (University of Victoria: Victoria, B.C.). Fall 2022.